


INTRODUCTION

The Seattle Laminating Adapter is designed for use exclusively in lower limb prostheses. This adapter can be easily incorporated into a two-part socket lamination for AK limbs. The section below indicates the maximum patient weight.

Product Code	Description	Patient Weight Limit
AAASS237	Stainless Steel, Anchor, Rotatable	136 kg / 300 lb (Med. Activity)
SCA237-01	3-Prong Laminating Dummy	N/A



LOW	Walking with aid
MED LOW	Limited walking
MED	Walking
MED HIGH	Jogging, light sports
HIGH	Running, basketball, farming, and other strenuous activities

TABLE 1: Activity Levels

CHART 1: AK lamination adapter weight limit selection chart.

INSTALLATION AND USE

Note: Recommended installation and use procedures must be followed for maximum safety and service life.

1. Begin lamination using 2 layers of Nyglass stockinette, 1 layer of carbon cloth, and 1 layer of fiberglass matting.
2. Apply 2 strips of carbon cloth — 1 from the medial to lateral areas of the socket and 1 wrapped around the socket just below the patella.
3. Apply 2 layers of Nyglass stockinette.
4. Continue laminating in the usual manor. After the lamination has cured, roughen out surface and align laminating adapter accordingly to the distal end.
5. Attach the adapter to the lamination using an acrylic bonding paste.
Note: Be sure not to leave any gaps or air pockets between the lamination and the adapter.
6. If applicable, place the pyramid cap onto the unit.
7. Proceed with final layup over the adapter using 1 layer of carbon cloth and 2 layers of fiberglass stockinette. The second layer of fiberglass should extend half the length of the socket. Finally, apply 2 layers of Nyglass stockinette.
8. Continue with lamination of remainder of socket using resin system of choice.
9. Fill the receptacle with putty or Trulife’s 3-Prong Laminating Dummy (SCA237-01). Remove the tightening screw and apply Vaseline or other grease before replacing.

10. Proceed with I-Beam layup over the adapter and continue with lamination of remainder of socket as normal.
11. Once lamination is complete, remove the putty or Trulife's 3-Prong Laminating Dummy. Clear the threads of any putty, grease, and/or resin.
12. Screw in the attachment unit until the attachment piece makes contact with the top of the flange unit, again ensuring that there is no foreign material between assemblies. The attachment piece may now be rotated in the reverse direction in order to obtain proper alignment. However, it is extremely important that the maximum amount of threads be engaged for full component strength.
13. Tighten the clamp ring attachment screw to 9 Nm (6.3 ft lb) and the set screws (if applicable) to 15 Nm (11 ft-lb). To maintain bolt tightness, apply Loctite 242 removable thread locking compound to the threads of the clamp bolt. Loctite requires several hours to cure completely.

Warning: Failure to follow the installation and use procedures set forth above may lead to structural failure of the components subjecting the user to a risk of serious personal injury.

QUESTIONS

Contact Customer Service in the U.S. at 888-878-1238, or fax 888-878-1237.

If calling from outside the U.S., contact Customer Service at 360-697-5656, or fax 360-697-6843.

Visit Trulife online at www.trulife.com.

Trulife has appointed Medical Device Safety Service (MDSS) of Hannover, Germany to act as our EU authorized representative. They may be contacted at:

MDSS GmbH
Schiffgraben 41
30175 Hannover
Germany
Phone (+49)-511-6262 8630
FAX (+49) -511-6262 8633

LIMITED WARRANTY

Trulife warrants that the AK Laminating Adapter will be free from defects in material and workmanship for two (2) years from the date of purchase.

This warranty will not apply if the product has been damaged by misuse, abuse, neglect, improper care, failure to follow instructions, abnormal wear and tear, or in the event that the adapter has been modified/repared by persons unauthorized by Trulife.

If a defect in material or workmanship is found during the warranty period, Trulife will, at Trulife's option, either repair or replace the product. If it is not possible to repair or replace the product, Trulife will be limited to refunding the purchase price.

Trulife will not be liable under any legal theory for any direct, indirect, special, incidental or consequential damages arising from the use of or inability to use this product

The application guidelines for this Trulife product are for the use of and by certified, qualified practitioner only. Patients are not to attempt to apply or adjust the item unless expressly instructed to do so by the practitioner responsible for the prescription and/or initial fitting of the device. All patient questions should be referred to the practitioner and not to the manufacturer. The manufacturer warrants only that the enclosed product has been inspected for quality and can be effective for certain indications, but final decisions and ongoing monitoring must be made by the prosthetic professional(s) prescribing and/or fitting the device to determine its effectiveness for an individual patient. Patient compliance is an integral part of the entire protocol and must be adhered to in order to avoid potential problems and to maximize the effectiveness of the prescribed product.

As a condition of the sale of any Trulife product, this prosthesis is restricted to a "Single Patient Use Only" by the originally fitted patient in order to protect the care provider and the patient against potentially adverse consequences of infectious disease transmission, material instability in adapting to the configuration of the original user and/or decrease in effectivity. Any express or implied warranties are voided if the product is reused or fitted to another patient. Additionally, a license of right to use under any relevant patents pertaining to the product is terminated with the cessation of use by the original patient. As with all Trulife products, this product must be prescribed and applied by a qualified practitioner to determine it meets the needs of the particular patient and accomplishes the desired results.