

TruCare

Pressurecare you can wear



Following extensive research and testing, **Trulife** have developed a new **TruCare** range which is unlike anything currently available in the field of pressurecare.

Designed to allow multifunction use both on the body and on materials that come in contact with pressure points. We believe these new products can enhance the care of patients at risk of suffering from pressure sores and ensure optimum comfort.

The range is designed to protect the skin from pressure. The combination of skin compatible silicone gel and adhesive provide the patient with excellent levels of protection and comfort.

There are currently four shapes available in this line.

TruCare



TruCare - Pressurecare you can wear
Ideal for Long-term / Elderly Care

Please note all products are sold in Pairs

TruCare

TC-01



The **Elbow** product can be applied directly onto the patients elbow and reapplied numerous times without losing tackiness. May also be useful for wheelchair users who may find the arm rest uncomfortable.



TC-02



The **Wrist** product can be applied directly onto the patient's wrist protecting the ulna and radius. May also be useful for wheelchair users who may find the arm rest uncomfortable.



TC-04



The **Positioning Pad**, is especially useful for persons in a wheelchair who may have muscle atrophy and an inclination to leaning on one side of the wheelchair. It can be applied directly onto the patient's body and reapplied numerous times. Can either be applied to the patient or a hard surface that the patient may come in contact with.



TC-03



The **Heel** product can be applied directly onto the heel and reapplied numerous times without losing tackiness, especially designed to redistribute the pressure on the calcaneus and the bottom of the Tibia and Fibula.



Applying Products

Elbow Application
Product Code: TC-01



Wrist Application
Product Code: TC-02



Heel Application
Product Code: TC-03



Positioning Pad Application
Product Code: TC-04



Instructions for General Use

Remove the pad from its protective backing film



Apply directly to skin



What Condition should skin be in before use?

Skin Should be clean and dry

Can the product be used on damaged skin?

The product should not be used on damaged skin

How long can a patient wear the product for?

A patient can wear the product for up to 5 hours

The product should be checked intermittently to ensure it has not become detached from the patient

Can the product be reused?

The product is for single person use

It can be reused multiple times provided it is cleaned between uses

How long can the product be used for?

Discontinue use when the pad loses its adhesion

Cleaning instructions

Make sure the product is clean before use.
Wash the pad with hot soapy water and dab dry with paper towel
before re-application.
Cleaning reinvigorates tackiness.

Storage

Do not store in direct sunlight.
Pad should be stored in its original packaging.



TruCare

Pressurecare you can wear

The bottom right corner of the image features decorative graphic elements consisting of several overlapping, curved lines. These lines are light gray and white, creating a sense of motion and depth. The lines curve from the bottom right towards the center of the page.