



## Trulife Lower Limb Systems Components

| ABOVE KNEE   | ZENITH  | ZENITH LP                  | TRIUMPH                    | ENERGY  | LIGHTFOOT <sup>2</sup>     | EDGE                           | KINETIC LIGHT                  | KINETIC                        | NATURAL                  | SACH                     | CHILD'S PLAY         |                      |
|--------------|---|----------------------------|----------------------------|---|----------------------------|--------------------------------|--------------------------------|--------------------------------|--------------------------|--------------------------|----------------------|----------------------|
|              | Impact Level                                      | K3/K4                      |                            |   | K3                         |                                | K2                             |                                |                          |                          | K1                   | K3/K4                |
| Weight Limit | 100 kg/220 lb (SSK615)<br>125 kg/275 (SSK615-THR) |                            |                            | 100 kg/220 lb (SSK615)<br>125 kg/275 (SSK615-THR) |                            | 100 kg/220 lb                  |                                |                                |                          | 102 kg/225 lb            | 55 kg/121 lb         |                      |
|              | AAASMP  | AAASMP                     | AAASMP                     | AAASMP  | AAASMP                     | AAASMP                         | AAASMP                         | AAASMP                         | AAASMP                   | AAASMP                   | AAAKSMP              |                      |
|              | AAA212  | AAA212                     | AAA212                     | AAA212  | AAA212                     | AAA212                         | AAA212                         | AAA212                         | AAA212                   | AAA212                   | AAAK211              |                      |
|              | AAA214-01<br>(SSK615 only)                        | AAA214-01<br>(SSK615 only) | AAA214-01<br>(SSK615 only) | AAA214-01<br>(SSK615 only)                        | AAA214-01<br>(SSK615 only) | AAA214-01                      | AAA214-01                      | AAA214-01                      | AAA214-01                | AAA214-01                |                      |                      |
|              | SSK615<br>or<br>SSK615-THR                        | SSK615<br>or<br>SSK615-THR | SSK615<br>or<br>SSK615-THR | SSK615<br>or<br>SSK615-THR                        | SSK615<br>or<br>SSK615-THR |                                |                                |                                |                          |                          |                      |                      |
|              |   |                            |                            |   |                            | SSK617                         | SSK617                         | SSK617                         |                          |                          |                      |                      |
|              |   |                            |                            |   |                            |                                |                                |                                | SSK602L                  | SSK602L                  |                      |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          | SSK610A              |                      |
|              | AAA200  | AAA200                     | AAA200                     | AAA200  | AAA200                     | AAA200<br>(Included with foot) | AAA200<br>(Included with foot) |                                | AAA200                   | AAA200                   |                      |                      |
|              |   |                            |                            |   |                            |                                |                                | AAA202<br>(Included with foot) |                          |                          |                      |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          | AAAK135              |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          | AAAK285              |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          | SFS24x<br>(Optional) |                      |
| Foot         | SZN200/203  | SZN100/103                 | STF190/193                 | SEF190/193  | SLF195/198                 | SKF300/303                     | SKF100/103                     | SKF200/203                     | SNF160/163<br>SNF170/173 | SCH110/113<br>SCH120/123 | SEF132/144           |                      |
| BELOW KNEE   | ZENITH  | ZENITH LP                  | TRIUMPH                    | ENERGY  | LIGHTFOOT <sup>2</sup>     | EDGE                           | KINETIC LIGHT                  | KINETIC                        | NATURAL                  | SACH                     | ZUMO                 | CHILD'S PLAY         |
|              | Impact Level                                      | K3/K4                      |                            | K3  |                            | K2                             |                                |                                | K1                       | K3                       | K3/K4                |                      |
| Weight Limit | 136 kg/300 lb                                     |                            | 160 kg/350 lb              | 136 kg/300 lb                                     |                            | 160 kg/350 lb                  |                                |                                | 102 kg/225 lb            | 205 kg/450 lb            | 65 kg/144 lb         |                      |
|              | AAASMP  | AAASMP                     | AAASMP                     | AAASMP  | AAASMP                     | AAASMP                         | AAASMP                         | AAASMP                         | AAASMP                   | AAASMP                   | AAASMP               | AAAKSMP              |
|              | AAA210  | AAA210                     | AAA210                     | AAA210  | AAA210                     | AAA210                         | AAA210                         | AAA210                         | AAA210                   | AAA210                   | AAA210               | AAAK210              |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          | AAAHDTi210           |                      |
|              | AAA225  | AAA225                     | AAA225                     | AAA225  | AAA225                     | AAA225                         | AAA225                         | AAA225                         | AAA225                   | AAA225                   | AAA225               |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          | AAAHDTi225           |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          |                      | AAAK100              |
|              | AAA200<br>or<br>AAA135                            | AAA200<br>or<br>AAA135     | AAA200<br>or<br>AAA135     | AAA200<br>or<br>AAA135                            | AAA200<br>or<br>AAA135     | AAA200<br>(included with foot) | AAA200<br>(included with foot) |                                | AAA200<br>or<br>AAA135   | AAA200<br>or<br>AAA135   |                      |                      |
|              |   |                            |                            |   |                            |                                |                                | AAA202<br>(included with foot) |                          |                          |                      |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          | AAAHDTi200           |                      |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          |                      | AAAK135              |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          |                      | AAAK285              |
|              |   |                            |                            |   |                            |                                |                                |                                |                          |                          |                      | SFS24x<br>(Optional) |
| Foot         | SZN200/203  | SZN100/103                 | STF190/193                 | SEF190/193  | SLF195/198                 | SKF300/303                     | SKF100/103                     | SKF200/203                     | SNF160/163<br>SNF170/173 | SCH110/113<br>SCH120/123 | SZM450/453           | SEF132/144           |



# Trulife AAA Modular Components

| A A A A L U M I N U M A D A P T E R S             |           |  |                        |                    |                           |                |               |  |
|---|-----------|--|------------------------|--------------------|---------------------------|----------------|---------------|--|
| IMAGE   | MODEL NO. | DESCRIPTION                                  | MATERIAL               | CONNECTION         | PRODUCT DIMENSIONS        | PRODUCT WEIGHT | WEIGHT LIMIT  | INCLUDES   |
|   | AAASS237  | 3-Prong Rotatable Laminating Piece           | Stainless Steel        | M36 Thread         |                           | 107 g/3.8 oz   | 136 kg/300 lb | Clamp Screw  |
|   | AAA210    | 4-Hole Pyramid Adapter                       | Aluminum               | Male, Four-Hole    | 10 mm ID                  | 90 g/3.2 oz    | 160 kg/350 lb | Mounting Screws                                      |
|   | AAA211    | 4-Hole Female Adapter                        | Aluminum               | Female, Four-Hole  |                           | 88 g/3.1 oz    | 160 kg/350 lb | Set Screws<br>Mounting Screws                        |
|   | AAA212    | 4-Hole Rotatable Flange                      | Aluminum               | Four-Hole          |                           | 58 g/2.0 oz    | 160 kg/350 lb | Mounting Screws                                      |
|   | AAA214-01 | Rotatable Female Insert                      | Aluminum               | Female, M36 Thread |                           | 67 g/2.3 oz    | 160 kg/350 lb | Set Screws   |
|   | AAA216-01 | Rotatable Male Pyramid Insert                | Aluminum               | Male, M36 Thread   | 10 mm ID                  | 54 g/1.9 oz    | 160 kg/350 lb |  |
|   | AAA208    | 4-Hole Double Offset Pyramid Adapter         | Aluminum               | Male, Four-Hole    | 10 mm ID                  | 90 g/3.2 oz    | 160 kg/350 lb | Mounting Screws                                      |
|   | AAA209    | 4-Hole Single Offset Pyramid Adapter         | Aluminum               | Male, Four-Hole    | 10 mm ID                  | 90 g/3.2 oz    | 160 kg/350 lb | Mounting Screws                                      |
|   | AAA100    | Adjustable Height Clamp Adapter              | Aluminum               | Female             | 30 mm ID<br>75 mm Height  | 120 g/4.2 oz   | 160 kg/350 lb | Set Screws<br>Clamp Screw                            |
|   | AAA225    | Clamp Adapter                                | Aluminum               | Female             | 30 mm ID<br>48 mm Height  | 90 g/3.2 oz    | 160 kg/350 lb | Set Screws<br>Clamp Screw                            |
|   | AAA221    | Double-Ended Adapter                         | Aluminum               | Female             | 35 mm Height              | 109 g/3.8 oz   | 160 kg/350 lb | Set Screws   |
|   | AAA222    | Double-Ended Adapter                         | Aluminum               | Female             | 45 mm Height              | 125 g/4.4 oz   | 160 kg/350 lb | Set Screws   |
|   | AAA223    | Double-Ended Adapter                         | Aluminum               | Female             | 60 mm Height              | 142 g/5.0 oz   | 160 kg/350 lb | Set Screws   |
|   | AAA224    | Double-Ended Adapter                         | Aluminum               | Female             | 75 mm Height              | 158 g/5.6 oz   | 160 kg/350 lb | Set Screws   |
|   | AAA285    | Foot Pyramid                                 | Aluminum               | Male, M10 Thread   | 24 mm Height<br>51 mm A-P | 67 g/2.4 oz    | 136 kg/300 lb |  |
|   | AAASMP    | Mounting Plate                               | Aluminum               | Four-Hole          | 17 mm Height              | 118 g/4.2 oz   | 205 kg/450 lb | Mounting Screws<br>Expulsion Valve<br>Fab Components |
|   | AAA200    | Aluminum Pylon Bonded to Female Adapter      | Aluminum               | Female             | 30 mm<br>420 mm Length    | 324 g/11.4 oz  | 160 kg/350 lb | Set Screws   |
|   | AAA202    | Aluminum Pylon Bonded to Male Adapter        | Aluminum               | Male               | 30 mm<br>420 mm Length    | 329 g/11.6 oz  | 160 kg/350 lb |  |
|   | AAA135    | Variable Height Adapter                      | Aluminum               | Female             | 30 mm<br>178 mm Length    | 190 g/6.7 oz   | 160 kg/350 lb | Set Screws   |
|   | AAA201    | Pylon  | Aluminum               |                    | 30 mm<br>428 mm Length    | 308 g/10.9 oz  | 160 kg/350 lb |  |
| A A A H D T i H E A V Y - D U T Y A D A P T E R S |           |  |                        |                    |                           |                |               |  |
| IMAGE   | MODEL NO. | DESCRIPTION                                  | MATERIAL               | CONNECTION         | PRODUCT DIMENSIONS        | PRODUCT WEIGHT | WEIGHT LIMIT  | INCLUDES   |
|   | AAAHTi210 | 4-Hole Pyramid Adapter                       | Titanium               | Male, Four-Hole    |                           | 110 g/3.5 oz   | 205 kg/450 lb | Mounting Screws                                      |
|   | AAAHTi217 | Splined 4-Hole Rotatable Adapter             | Titanium               | Male, Four-Hole    |                           | 128 g/4.5 oz   | 205 kg/450 lb | Mounting Screws                                      |
|   | AAAHTi218 | Splined 4-Hole Rotatable Adapter             | Titanium               | Female, Four-Hole  |                           | 138 g/4.9 oz   | 205 kg/450 lb | Set Screws<br>Mounting Screws                        |
|   | AAAHTi225 | Adjustable Height Clamp Adapter              | Titanium               | Female             | 30 mm ID<br>50 mm Height  | 102 g/3.6 oz   | 205 kg/450 lb | Set Screws<br>Clamp Screw                            |
|   | AAAHTi221 | Double-Ended Adapter                         | Titanium               | Female             | 35 mm Height              | 118 g/4.2 oz   | 205 kg/450 lb | Set Screws   |
|   | AAAHTi200 | Bonded Pylon, 16.5"                          | Titanium/<br>Composite | Female             | 30 mm OD<br>42 cm Length  | 328 g/11.6 oz  | 205 kg/450 lb | Set Screws   |
| A A A K C H I L D ' S P L A Y A D A P T E R S     |           |  |                        |                    |                           |                |               |  |
| IMAGE   | MODEL NO. | DESCRIPTION                                  | MATERIAL               | CONNECTION         | PRODUCT DIMENSIONS        | PRODUCT WEIGHT | WEIGHT LIMIT  | INCLUDES   |
|   | AAAK210   | Child's Play 4-Hole Pyramid Adapter          | Aluminum               | Male, Four-Hole    |                           | 49 g/1.7 oz    | 65 kg/144 lb  | Mounting Screws                                      |
|   | AAAK211   | Child's Play 4-Hole Female Adapter           | Aluminum               | Female, Four-Hole  |                           | 39 g/1.4 oz    | 65 kg/144 lb  | Set Screws<br>Mounting Screws                        |
|   | AAAK221   | Child's Play Double-Ended Adapter            | Aluminum               | Female             | 22 mm ID<br>32 mm Height  | 50 g/1.8 oz    | 65 kg/144 lb  | Set Screws   |
|   | AAAK100   | Child's Play Adjustable Height Clamp Adapter | Aluminum               | Female             | 22 mm ID                  | 58 g/2.0 oz    | 65 kg/144 lb  | Set Screws<br>Clamp Screw                            |
|   | AAAK225   | Child's Play Clamp Adapter                   | Aluminum               | Female             | 22 mm ID                  | 40 g/1.4 oz    | 65 kg/144 lb  | Set Screws<br>Clamp Screw                            |
|   | AAAK135   | Child's Play Variable Height Adapter         | Aluminum               | Female             | 22 mm<br>178 mm Length    | 87 g/3.1 oz    | 65 kg/144 lb  | Set Screws   |
|   | AAAK285   | Child's Play Energy Foot Pyramid             | Aluminum               | Male, M8 Thread    | 37 mm A-P<br>14 mm Height | 33 g/1.2 oz    | 65 kg/144 lb  |  |
|   | AAAKSMP   | Child's Play Mounting Plate                  | Aluminum               | Four-Hole          | 5.8 cm OD                 | 56 g/2.0 oz    | 65 kg/144 lb  | Mounting Screws<br>Fab Components                    |